
COVID-19 Ag FIA

STANDARD™ F COVID-19 Ag FIA

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

KIT CONTENTS


MATERIALS REQUIRED BUT NOT PROVIDED

- STANDARD F Analyzer
- Timer

SPECIMEN COLLECTION AND STORAGE

■ Nasopharyngeal swab


- To collect a nasopharyngeal swab specimen, insert the sterile swab into the nostril that present the most secretion under visual inspection.
- Using gentle rotation, push the sterile swab until resistance is met at the level of the turbinate.
- Rotate the sterile swab a few times against the nasopharyngeal wall.
- Remove the sterile swab from the nostril carefully.
- Specimen should be tested as soon as possible after collection.
- Specimen may be stored at room temperature for up to 1 hours or at 2-8°C/ 36-46°F for up to 4 hours in a clean, dry, closed container prior to testing.

■ Viral transport medium

The following viral transport media listed in Table 1 were validated by SD BIOSENSOR'S R&D Department and determined to be compatible with the STANDARD F COVID-19 Ag FIA. The viral transport media are divided into two different groups by manufacturing brand. Find a brand of viral transport medium you use, and select a correct group button while sample selecting phase.

Table. 1 Recommended Viral Transport Media

Group	Viral Transport Medium
VTM – Group 1	COPAN UTM™ Starswab® Multitrans™ System ASAN Transport Medium STANDARD™ Transport Medium BD Universal Viral Transport
VTM – Group 2	Noble BIO REST™ UTM FA Transport Medium


Only viral transport medium verified by SD BIOSENSOR can be used as specimen for STANDARD F COVID-19 Ag FIA. Please check the brand of viral transport medium and select a correct group button while sample selecting phase.

PREPARATION AND TEST PROCEDURE

■ Preparation

- Allow test device and collected specimen to room temperature (15-32°C / 59-86°F) prior to testing.
- Carefully read instructions for using the STANDARD F COVID-19 Ag FIA.
*Check the expiry date at the back of the foil pouch. Do not use the kit, if expiry date has passed.
- Check the test device and the desiccant.

PERFORMANCE CHARACTERISTICS

■ Clinical Evaluation


The sensitivity and specificity of the STANDARD F COVID-19 Ag FIA for qualitative detection of SARS-CoV-2 antigen was established in retrospective, single institute, randomized, single-blinded studies conducted at a trial sites in Republic of Korea during the 2020 SARS-CoV-2 pandemic. A total of 155 specimens were tested using the STANDARD F COVID-19 Ag FIA with STANDARD Analyzer. Specimens consisted of nasopharyngeal swabs from symptomatic patients. The performance of STANDARD F COVID-19 Ag FIA was compared to a commercialized molecular assay. The STANDARD F COVID-19 Ag FIA showed 89.09% of sensitivity and 96% of specificity. Summary of the sensitivity and specificity of STANDARD F COVID-19 Ag FIA compared to PCR.

	PCR		Total
	Positive	Negative	
STANDARD F COVID-19 Ag FIA	49	4	53
	Negative	96	102
	Total	55	100
Sensitivity	89.09% (49/55), (95% CI: 77.75%–95.89%)		
Specificity	96% (96/100), (95% CI: 90.07%–98.90%)		

ANALYTICAL PERFORMANCE


■ Limit of Detection (LOD): The study used "SARS-CoV-2 (2019-nCoV) NCCP 43326/2020 / Korea" strain. The titer of cultured virus was confirmed by PCR. The cell is inactivated and slipped into nasopharyngeal swab specimen. The LO is 7.81 X 10⁻³ TCID₅₀/ml.

2019-nCoV Strain Tested		NCCP 43326/2020 / Korea									
Stock 2019-nCoV Titer		1 X 10 ³ TCID ₅₀ /ml									
Serial 2-Fold Dilution		1/100	1/200	1/400	1/800	1/1600	1/3200	1/12800	1/25600	1/51200	
Concentration in Dilution Tested		1 X 10 ³	5 X 10 ²	2.5 X 10 ²	1.25 X 10 ²	6.25 X 10 ¹	3.13 X 10 ¹	7.81 X 10 ⁰	3.91 X 10 ⁰	1.95 X 10 ⁰	
Replicate 1	POS	POS	POS	POS	POS	POS	POS	POS	POS	NEG	NEG
Replicate 2	POS	POS	POS	POS	POS	POS	POS	POS	POS	NEG	NEG
Replicate 3	POS	POS	POS	POS	POS	POS	POS	POS	POS	NEG	NEG
Replicate 4	POS	POS	POS	POS	POS	POS	POS	POS	POS	NEG	NEG
Replicate 5	POS	POS	POS	POS	POS	POS	POS	POS	POS	NEG	NEG
Call Rate		100%	100%	100%	100%	100%	100%	100%	100%	20%	0%
Lowest Concentration with Uniform Positivity per Analyte		7.81 X 10 ⁻³ TCID ₅₀ /ml									
Limit of Detection (LOD) per Virus Strain		7.81 X 10 ⁻³ TCID ₅₀ /ml									

Do not write on the barcode or damage the barcode of the test device.


■ Extraction of specimen [Nasopharyngeal swab]

- Insert nasopharyngeal swab into an extraction buffer tube. While squeezing the buffer tube, stir the swab more than 5 times.
- Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.
- Press the nozzle cap tightly onto the tube.

[Specimen in viral transport medium]

- Collect 350µl of viral transport medium specimen using a micropipette.

■ Analysis of specimen

STANDARD TEST mode STANDARD F100, F200 and F2400 analyzer

- Take the test device out of the foil pouch and place it on a flat and dry surface. Write patient information on the label of test device.

STANDARD F2400 analyzer

Workplace -> *Run Test* -> Insert patient ID and / or operator ID on the analyzer

STANDARD F100 and F200 analyzer


Standard Test mode -> Insert patient ID and / or operator ID on the analyzer

- Insert the test device to the test slot of the analyzer. When inserting the test device to the analyzer, the analyzer will read the barcode data, and check the test device is valid.

- Select correct specimen type (Nasopharyngeal swab / VTM – Group 1 / VTM – Group 2). When using specimen in the viral transport medium, a correct group button should be selected for adequate test result.

- Apply 4 drops of extracted specimen to the specimen well in the test device.

STANDARD F100 and F200 analyzer

READ ONLY mode

- Take the test device out of foil pouch and place it on a flat and dry surface. Write patient information on the label of test device.

- Apply 4 drops of extracted specimen to the specimen well in the test device.

- Incubate the test device for 30 minutes outside of the analyzer.

- Prepare a STANDARD F Analyzer and select the 'Read Only' mode according to the analyzer's manual.

- Insert the test device to the analyzer.

- Select correct specimen type (Nasopharyngeal swab / VTM – Group 1 / VTM – Group 2). When using specimen in the viral transport medium, a correct group button should be selected for adequate test result.

- The analyzer will automatically scan and display the test result immediately after specimen type selection.

Positive results should be considered in conjunction with the clinical history and other data available to the physician.

INTERPRETATION OF TEST RESULT

Result	COI (Cutoff Index) value	SARS-COV-2 Ag
Positive	COI ≥ 1.0	Positive for SARS-COV-2 Ag
Negative	COI < 1.0	Negative for SARS-COV-2 Ag
Invalid	COI value is not displayed	Retest should be performed with a new test device and a new patient's specimen.

The test result of a specimen is given either as Positive (+) / Pos (+) or Negative (-) / Neg (-) with a COI (cut off index) value. The COI is a numerical representation of the measured fluorescence signal.

EXPLANATION AND SUMMARY

■ Introduction

Coronavirus is a single-stranded positive-sense RNA virus with an envelope of about 80 to 120 nm in diameter. Its genetic material is the largest of all RNA viruses and is an important pathogen of many domestic animals, pets, and human diseases. It can cause a variety of acute and chronic diseases. Common signs of a person infected with a coronavirus include respiratory symptoms, fever, cough, shortness of breath, and dyspnea. In more severe cases, infection can cause pneumonia, severe acute respiratory syndrome, kidney failure, and even death. The 2019 new coronavirus, or '2019-nCoV', was discovered because of Wuhan Viral Pneumonia cases in 2019, and was named by the World Health Organization on January 12, 2020, confirming that it can cause colds and the Middle East Respiratory Syndrome (MERS) and more serious diseases such as acute respiratory syndrome (SARS). This kit is helpful for the auxiliary diagnosis of coronavirus infection. The test results are for clinical reference only and cannot be used as a basis for confirming or excluding cases alone.

■ Intended use

STANDARD F COVID-19 Ag FIA is the fluorescent immunoassay for the qualitative detection of specific nucleoprotein antigens to SARS-CoV-2 present in human nasopharynx. STANDARD F COVID-19 Ag FIA should be used with the STANDARD F Analyzers manufactured by SD BIOSENSOR. This test is for in vitro professional diagnostic use and intended as an aid to early diagnosis of SARS-CoV-2 infection in patient with clinical symptoms with SARS-CoV-2 infection. It provides only an initial screening test result. More specific alternative diagnostic methods should be performed in order to obtain the confirmation of SARS-CoV-2 infection.

■ Test principle

STANDARD F COVID-19 Ag FIA is based on immunofluorescence technology with STANDARD F Analyzer to detect SARS-CoV-2 nucleoproteins. STANDARD F COVID-19 Ag FIA has a test line which is coated with monoclonal anti-SARS-CoV-2 antibody. The patient's specimen is applied into the specimen well along the test device and the specimen migrates through the membrane. If SARS-CoV-2 viral antigen is present in the patient's specimen, it will react with europium conjugated monoclonal anti-SARS-CoV-2 antibody in the conjugation pad and form antibody-antigen fluorescence particle complex. This complex move along to the membrane to be captured by the anti-SARS-CoV-2 antibody on the test line and make fluorescence signal. The intensity of the fluorescence light generated on the membrane is scanned by the STANDARD F Analyzer manufactured by SD BIOSENSOR. STANDARD F Analyzer can analyze the presence of the SARS-CoV-2 in the clinical specimen by processing the results using pre-programmed algorithms and display the test result on the screen.

KIT STORAGE AND STABILITY

Store the kit at 2-30°C / 36-86°F out of direct sunlight. Kit materials are stable until the expiration date printed on the outer box. Do not freeze the kit.

WARNINGS AND PRECAUTIONS

- Do not re-use the test kit.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly after the tests are done.
- Clear up spills thoroughly using an appropriate disinfectant.
- Handle all specimens as if they contain infectious agents.
- Observe established precautions against microbiological hazards throughout testing processes.
- Dispose of all specimens and materials used to perform the test as bio-hazard waste. Laboratory chemical and biohazard wastes must be handled and discarded in accordance with all local, state, and national regulations.
- Desiccant in foil pouch is to absorb moisture and keep humidity from affecting products. If the moisture indicating desiccant beads change from yellow to green, the test device in the pouch should be discarded.
- Immediately use the test device after taking out of aluminum foil pouch.
- As the detection reagent is a fluorescent compound, no visible results will form on the test device.
- The barcode of the test device is used by analyzer to identify the type of test being run and to identify the individual test device so as to prevent to a second read of the test device by the same analyzer.
- Once a test device has been successfully scanned by analyzer, do not attempt to scan the test device again in the same analyzer.
- Improper specimen collection, handling or transport may yield inaccurate results.
- Do not write on the barcode or damage the barcode of the test device.

LIMITATION OF TEST

- The test procedure, precautions and interpretation of results for this test must be followed strictly when testing.
- This test detects the presence of SARS-CoV-2 in the specimen and should not be used as the sole criteria for the diagnosis of SARS-CoV-2 infection. Test results must be considered with other clinical data available to the physician.
- For more accuracy of immune status, additional follow-up testing using other laboratory methods is recommended.
- Neither the quantitative value nor the rate of SARS-CoV-2 concentration can be determined by this qualitative test.
- Failure to follow the test procedure and interpretation of test results may adversely affect test performance and/or produce invalid results.
- A negative result may occur if the concentration of antigen or antibody in a specimen is below the detection limit of the test or if the specimen was collected or transported improperly, therefore a negative test result does not eliminate the possibility of SARS-CoV-2 infection, and should be confirmed by viral culture or an molecular assay or ELISA.
- Positive test results do not rule out co-infections with other pathogens.
- Negative test results are not intended to rule in other coronavirus infection except the SARS-CoV.
- Children tend to shed virus for longer periods of time than adults, which may result in differences in sensitivity between adults and children.

QUALITY CONTROL

■ STANDARD F Analyzers calibration check


The calibration set test of STANDARD F Analyzers should be conducted according to the analyzer's manual.

[When to use calibration set]

- Before using the analyzer for the first time
- When you drop the analyzer
- Whenever you do not agree with your result
- When you want to check the performance of an analyzer and test device

[How to use calibration set]

- The calibration set test is a required function that ensures optimal performance by checking the internal analyzer optics and functions.
- Select the 'Calibration' menu.
- The specific calibration set is included with the analyzer.
- Insert the CAL-1 first, and then insert the CAL-2 for UV-LED testing and the CAL-3 for RGB-LED testing in order.

■ Internal procedural control

- The internal procedural control zone is in the end of the membrane of the test device. STANDARD F Analyzers read the fluorescence signal of the internal procedural control zone and decide whether the result is valid or invalid.
- The invalid result denotes that the fluorescence signal is not within the pre-set range. If the screen of STANDARD F analyzers shows 'Invalid Device', turn off and turn on the analyzer again and re-test with a new test device.

■ External quality control

- Positive and negative controls may be supplied with each kits or can be purchased from the distributors.
- It is recommended that positive and negative controls be run:
 - once for each new lot,
 - once for each untrained operator,
 - as required by test procedures in this instructions and in accordance with local, state and federal regulations or accreditation requirements.

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